

Image-Pro License Server System Requirements

Hardware Requirements

- ✓ Disk—500 MB
- ✓ RAM—4 GB
- ✓ CPU—2 GHz - 2 Cores

Browser Support

- ✓ Mozilla Firefox v43 & later
- ✓ Google Chrome v47 & later
- ✓ Microsoft Internet Explorer v10 & 11
- ✓ Microsoft Edge v38.1

Supported Platforms (CPU x86-64)

- ✓ Linux CentOS 7.0*
- ✓ Linux CentOS 8.0*
- ✓ Linux Ubuntu 20.04 LTS*
- ✓ Windows 8.1 and 10**
- ✓ Windows Server 2012 & 2012 R2**
- ✓ Windows Server 2016**
- ✓ Windows Server 2019**

* The FlexNet Embedded license server has only been validated on CentOS 7 and 8, but can also be run on RedHat Enterprise Linux 7 or 8.

** .NET Framework 4.5 or later is required for Windows platforms.

Virtual Machine (VM) Support

- ✓ VMware ESXi 6.5
- ✓ VMware Workstation 14.1.5
- ✓ Microsoft Hyper-V 6.3 on Windows Server 2016
- ✓ Citrix XenServer 7.2
- ✓ Oracle VirtualBox 6.1.18
- ✓ QEMU-KVM 2.7
- ✓ Parallels 15.1.2
- ✓ Google Compute Cloud
- ✓ Amazon EC2

Java

For machine where the FlexNet Embedded local license server is installed:

- ✓ Oracle Java SE 8, OpenJDK 8, or OpenJDK 11
- ✓ A 64-bit JRE for a 64-bit license server or a 32-bit JRE for a 32-bit license server (not adhering to this requirement can cause the license server to fail to start)